


VISICOIL™



Flexible soft tissue Marker for IGRT, SBRT & PT


**See it !
Trust it !
Treat it !**



See it !

VISICOIL™

Improved Marker Visibility by Design




1) Visibility during Marker Placement

- ✓ **More accurate marker placement:**
With TRUS / EUS / EBUS through highly echogenic VISICOIL surface


2) Better Visibility for Better Target Definition

- ✓ **VISICOIL aids in accurate target delineation particularly when tightening margins.**
CT: Reduced artifact for improved contouring
MRI: Compatible and visible with MRI for planning & fusion ¹⁾




3) Better Visibility for More Precise Image-Guided Patient Setup

- ✓ **Linear profile:**
Non-anatomic structure provides easy detection
- ✓ **Flexible marker design:**
Conforms with tissue motion




VISICOIL is Compatible with all Common Imaging Modalities




Ultrasound – Prostate




CBCT – Prostate




X-ray – Lung



MRI – Prostate



PET/CT – Pancreas




MV – Prostate

1) Jonsson JH, Garpebring A, Karlsson MG, et al. Internal fiducial markers and susceptibility effects in MRI-simulation and measurement of spatial accuracy. Int J Radiat Oncol Biol Phys 2012;1612-1618.

Trust it !

VISICOIL – Ideal for IGRT, SBRT, Proton Therapy and Robotic Radiosurgery



” The VISICOIL fiducials did not migrate and are easily identified on kilovoltage fluoroscopy using modern on-board imaging. ”

Creach KM, Azar R, Bradley JD, et al. Tumor localization with Visicoil ultra-thin fiducials for abdominal radiation therapy. Int J Radiat Oncol Biol Phys 2009;75:S589

” In the current study, fiducials implanted into a lung tumor showed no substantial migration during the entire course of treatment. ”

Kupelian PA, Forbes A, Willoughby TR, et al. Implantation and stability of metallic fiducials within pulmonary lesions. Int J Radiat Oncol Biol Phys 2007;777-785.

” The type of ...fiducial marker that we are employing has several advantages. One is the flexibility of the coiled marker. The other is the length. Together these characteristics allow ease of radial placement. We feel the coiled design lessens potential for migration. ”


McGuire FR, Kerley JM, Ochran T, et al. Radiotherapy monitoring device implantation into peripheral lung cancers: A therapeutic utility of electromagnetic navigational bronchoscopy. J Broncol 2007;189-192.

Treat it !

VISICOIL's Proven Efficiency and Ease of Use Enables you to Focus on your Patients

Faster, Easier and Safer Marker Placement

- ✓ Less needles = Less measuring = Faster implantations*
- ✓ Pre-loaded needles or loose VISICOIL = Easy to use
- ✓ VISICOIL products are stocked = Same day shipment available



*see implantation video at
YouTube (ibadosimetry)

Patient Setup Efficiency & Reliability

- ✓ Faster Setup through unambiguous marker detection
- ✓ VISICOIL is indicated for use to radiographically mark soft tissue for future therapeutic procedures

VISICOIL Support

- ✓ Dedicated team specializing in fiducial markers
- ✓ Local availability with vast experience
- ✓ Support your team through every step including placement and use
- ✓ As a courtesy, our partner Revenue Cycle is available for any reimbursement questions

EDUCATION

TRAINING

SUPPORT

NETWORKING



Innovations

THINNER NEEDLE
BETTER FOR PATIENT

Slim Line

PROVEN MARKER | SMALLER NEEDLE

The Proven 0.5 mm VISICOIL™ in Pre-loaded 19 Gauge Needle:

- ✓ Improved patient comfort during placement
- ✓ Designed to minimize edema
- ✓ Increased implantation safety especially for lung cases
- ✓ Recommended for: KV, CBCT, OBI, Ultrasound

= **30% Reduction in Needle Size, when
Compared to Typical Gold Seeds**



NEW!

Twin Line


PROVEN MARKER | TWO IN ONE

2 VISICOIL™ Markers with 1 Needle!

- ✓ Consistent marker spacing
- ✓ Half of the needles needed
- ✓ Faster and safer implantation



1 Twin-Line marker consisting of 2 pre-spaced
VISICOIL™ markers (0.5 cm each) with
radio-transparent spacer



2 Twin-Line VISICOIL™ markers
providing efficiently 4 reference points

VISICOIL™

What you see is what you treat!

VISICOIL is the original flexible linear marker for use with today's most accurate image guided radiotherapy treatment plans.



For more information or to order VISICOIL, please contact us:

Tel.: 901 432 7202

Toll Free: 866 429 0922

Fax: 901 432 7206

Email: VISICOIL@iba-group.com



www.VISICOIL.com

Protect,
enhance
and save
lives

iba